

A Closer Look at Compensation Methods in FlowJo

Amaya.Willemsen@cchmc.org



Compensation

- 4 different methods of compensation on FlowJo
 - Conventional
 - Autospill
 - “Spectral”
 - “Spectral” Autospill
- Compensation starts with good controls
 - Bright or brighter controls
 - Record plenty of events when acquiring
 - ~30,000 events are best
 - Or at least 100 events in the positive gate

What is Compensation





Figure 1. Accurate compensation aligns the MFI of the positive population to the MFI of the negative population.

- A correction of spillover of a fluorochrome into a secondary detector.
- Single Color Compensation controls are acquired and a positive and negative emission is defined



What is Autospill?

- A Flow Jo compensation algorithm that uses gate-less compensation
- Do not need to define a positive and negative population and do not need an unstained control
- Algorithm will define positive and negative populations
- Can use an unstained control to remove autofluorescence
 - Utilizing an empty detector that has the highest Autofluorescence
- Simple, and more efficient
- Helps remove biased positive and negative gating and human error
- [AutoSpill](#): a method for calculating spillover coefficients in high-parameter flow cytometry.
Authors: Carlos P. Roca, Oliver T. Burton, Teresa Prezzemolo, Carly E. Whyte, Richard Halpert, Łukasz Kreft, James Collier, Alexander Botzki, Josef Spidlen, Stéphanie Humblet-Baron, Adrian Liston.

What is Spectral Compensation?

- Need more detectors than fluorescent probes
- Utilizing spectral compensation, signal from non-primary detectors are still used to help unmix the true signal of the primary markers
- Same “process” as traditional compensation
 - i.e. define negative and positive population
- FlowJo spectral algorithm has an optimize weights function to utilize selected comp controls to compare the hypothetical spillover spread matrix (SSM) for different weights



Spillover Spread Matrix (SSM)

- A matrix that reveals how much fluorescence can spillover into each

SSM (Spillover Spread Matrix)

SSM	APCA	APC-Cy7-A	Alexa700	BB515-A	BB630-A	BB660-A	BB700-A	BB750-A	BB790-A	BUV396-A	BUV496-A	BUV563-A	BUV615-A	BUV661-A	BUV737-A	BUV805-A	BV421-A	BV480-A	BV570-A	BV605-A	BV650-A	BV711-A	BV750-A	BV786-A	PE-Cy5-A	PE-A	PE-Cy7-A	PECF504-A	
APC		1.12338	2.79735	0	0	1.3302	0.803615	0.670498	0.619803	0	0	0	0.178153	0.985026	0.857515	0.773112	0.22534	0.214154	0	0.156879	0.75195	0.390442	0.508582	0.473061	1.91928	0.494087	1.32326	0.394357	
APC-Cy7			4.84104	0.114748	0.275282	1.30188	0.354791	0.962837	1.26312	0	0	0	0.078791	0.272807	0.446186	1.40597	0.221141	0.209183	0	0.09877	0.217063	0.104088	0.198948	0.536302	0.581168	2.90521	0	0	
ALEXA 700	0.396896	7.24766		0.160912	0.780831	0.088211	1.28994	0.826428	0.701424	0.378361	0.326972	0.392132	0.240496	0.465403	1.13586	0.936619	0.367275	0.100224	0	0.324172	0.283231	0.367619	0.898561	0.767887	0.298943	0.356408	1.03066	1.69958	
APC-R700	0.754489	1.65489		0	0	0.257898	0.729263	0.662329	0.740682	0	0	0	0.127628	0.27043	0.664137	0.735164	0.121695	0	0.077977	0	0.190828	0.374801	0.508083	0.529383	0.354627	0.148525	1.05296	0.138185	
R728	0.35716	1.70005		0	0	0.178756	1.39334	0.468735	0.671072	0	0.114415	0.014451	0.054613	0	0.910907	1.03389	0.245453	0.157448	0	0.048675	0.055765	0.256048	0.292609	0.45253	0.21402	0.197117	1.36359	0.291147	
BB515	0.100933	0	0.031117		0.316454	0.190924	0.350283	0.221276	0.122298	0.06317	0	0	0.174565	0.141501	0.093383	0	0.049224	0.114584	0.084892	0	0.082863	0.076869	0.072632	0.085968	0.090007	0.129152	0	0.244821	
FITC	0.286658	0	0.286175	0.324821		0.29129	1.00081	0.43891	0	0.121852	0	0	0	0	0.202008	0.2466	0.148759	0	0.124917	0	0.004026	0.257063	0.064361	0.102587	1.18234	0.287193	0.564943	0	
BB630	1.15051	0.12168	0.789114	0.267795		1.32726	3.34131	1.58734	1.05972	0.303429	0.378902	1.31447	1.28345	0.577988	0.62196	0.290308	0.279527	0.356261	1.4503	1.11921	0.777884	0.791841	0.617019	0.852678	1.40413	0.556807	4.24614	0	
BB660	3.32895	0.907467	2.76473	0.357319	0.596821		4.75842	3.25496	1.89625	0.196445	0	0.094823	0.066271	1.51599	1.0884	1.14556	0.310154	0	0.247697	0.45597	2.04916	0.880264	1.11003	1.18979	1.34208	0.677105	0.410813	0	
BB700	1.01541	0.821174	1.94109	0.111335	0.240808	0.760508		2.82334	1.60795	0.073856	0.092748	0	0.142601	0.770888	0.967672	1.02485	0.283211	0.118864	0	0.244452	0.534204	1.14745	1.22746	0.926061	0.473054	0.215648	0.485194	0.217913	0
PercepcYS 5	1.18701	0.766202	1.77797	0	0.770639	0.921211		2.39418	1.46604	0.732747	0.383546	0	0.202836	0.841925	1.33352	0.789544	0.683996	0.656846	0	0.667995	1.09371	1.42162	1.0438	1.07717	0	1.69706	0.856803	0	
BB750	0.305512	0.900403	0.499555	0.105607	0.683832	0.486058	1.51678		2.4958	0.134812	0	0.068463	0.162616	0.194952	0.93806	1.27267	0.126372	0.176305	0	0.200059	0.226814	0.289772	0.934795	1.15906	0.078994	0.646534	0.658332	0.210068	0
BB790	0.397486	1.55059	0.431554	0.502883	2.93846	0.504944	1.96073	5.70179		0.040993	0	0	0.338315	0.169006	0.360766	0.523966	2.52391	0.229994	0.077289	0.426989	0.484158	0.335106	0.245335	0.496571	1.76624	0.298526	0.533629	0.64194	1.13004
BUV396	0.003308	0	0.202199	0	0	0	0.631855	0.203845	0.321495	0	0	0	0.425495	0.190813	0.421612	0	0.445388	0	0	0.242091	0	0.007064	0.291539	0.176598	0.446361	1.05749	0	1.02048	
BUV496	1.2885	0	0.792179	1.42407	1.45172	0	1.06526	0.026276	0	0	0	2.23673	2.07201	1.76499	0.794528	1.44681	0.59313	0.70011	0	0.831655	0.680161	0.709066	0.664638	0	0.675335	1.76268	0.041399	4.24479	
BUV563	1.22464	0	0.481878	0	3.10455	1.05133	1.22457	0.661454	0.377893	0	0	0	2.20005	1.88918	0.898665	0.789217	0.439116	0	0.396326	0.771334	0.423409	0.128017	0	0.100237	1.14336	6.53833	0.568341	5.91746	
BUV615	1.5681	0.320231	0.581181	0.031773	4.71004	1.07623	1.0812	0.736944	0.778387	0.058993	0	0	0.680092	2.63912	1.31237	1.34944	0.159915	0.193105	0.161917	1.37031	0.559905	0.335843	0.328682	0.322224	1.60685	3.40619	1.34454	10.2304	
BUV661	4.46763	1.01127	2.63628	0	0.693502	0.904923	0.773446	0.920997	0.594671	0	0.108803	0	0.406518	2.34698	1.88203	0.152262	0.203667	0.013008	0.228434	0.755777	0.415433	0.479853	0.441721	1.36599	0.662395	0.889524	0.956419	0	
BUV737	0.195628	1.53758	2.31896	0	0.181272	0.186895	1.48307	2.92234	1.21707	0.125723	0	0	0	0.196282	0	0	3.12146	0.241494	0.11185	0	0	0.109551	0.223371	0.603238	0.595978	0	0	0.614622	0.293702
BUV805	0	0	1.6189	0.136908	0	0	0.237752	0.314431	0.445214	0.287505	0	0	0.018437	0.067761	0.143057	0.250565	0	0	0.060965	0	0.056272	0.123254	0.223076	0.094884	0.315993	1.46095	0.050633	0	
BV421	0	0	0	0	0	0	0.284277	0	0	0.114093	0.286818	0	0	0.166907	0.008481	0.135456	0	0	0.493474	0	0.197802	0.275888	0.079887	0.08381	0	0.140682	0.596238	0.079083	0
BV480	0.153152	0	0.337346	0.727896	1.18873	0.300407	0.834699	0.169692	0	0	0.613169	0.864527	0.452831	0.274219	0.457752	0.376961	0	0	0.707186	0.987485	0.844785	0.35426	0.65046	0.422048	0.282155	1.12787	0	1.14091	
BV510	1.23495	0	1.55036	1.11024	3.1996	0.463732	1.0619	0.61631	0	1.1036	0	0	2.35529	2.98378	2.74552	1.87959	2.42546	1.88041	0	2.11873	3.22036	2.10138	1.47773	1.88419	1.84127	0.507465	5.18644	0.698125	4.57556
BV570	0.755939	0	1.39863	0.477006	3.60417	1.08665	1.78135	0.910348	0.743261	0.817578	0	0	1.87368	2.35654	1.19625	0.673571	0.769706	1.14305	0.465965	2.38559	1.63324	0.993579	1.21989	1.14934	1.60369	7.90481	1.19954	5.71187	
BV605	1.20622	0	1.61392	0.113303	2.96814	1.39945	1.8128	1.0269	0.898702	0	0.497792	0.644784	2.1114	1.97027	1.40232	1.37558	0.875495	0.36455	0.809735	1.73415	1.50174	1.49862	1.35594	1.65376	3.99091	1.4435	5.48146	0	
BV650	2.87813	0.670604	2.95202	0	1.24349	0.998393	1.37088	0.88789	0.58999	0	0	0.206695	0.795475	2.25869	1.97326	1.8872	0.633767	0.218025	0.210372	0.66637	0	1.73125	2.29725	1.76018	1.40043	0.535761	0.834351	1.34372	
BV711	0.634265	1.5459	6.37081	0	0	0.495586	2.5346	2.22672	1.4999	0.247658	0.149611	0.166731	0.261129	0.596021	3.08119	3.99235	0.907932	0	0	0.111704	0.632433	0	2.7216	2.81776	0.330793	0.444456	0.506933	0	
BV750	0.188165	1.70392	1.13566	0	0	0	0.926882	2.4507	1.06534	0.10133	0	0	0.1194374	0.134973	2.98437	3.2735	0.727731	0.416566	0.21411	0.171638	0.252636	0.508185	0	0.25169	0.865508	0.198956	0		
BV786	0.202576	2.30054	0.427883	0.538378	0.49541	0.061688	0.352157	0.823412	1.14645	0.4516	0.366356	0	0.013265	0.20964	1.04918	4.3041	2.5015	0.287371	0.282563	0.185761	0.224465	0.125932	0.922329	0	0.155739	0.082977	0.99322	0.256066	
PCy5	2.87539	0.790363	2.02519	0	0.270387	2.91966	0.78744	3.42573	1.98089	0.161082	0	0	1.20197	0.671505	0.585297	0.192019	0	0	0.2128	1.40165	0.623618	0.756443	0.74354	0	0.623937	1.95075	0.472786	0	
PE	0.506759	0	0.261808	0	1.69286	0.713827	1.49376	0.808046	0.49709	0.205288	0.053004	0.765737	0.630142	0.560914	0.29608	0.191519	0.428046	0.111375	0.904157	1.02724	0.497073	0.307193	0.351771	0.281128	0.707061	0	0.392039	2.44152	
RYG584	0.251521	0	0.276246	0.00611	1.74675	0.571513	0.953132	0.511598	0.395956	0.062302	0.054105	0.112575	0.148335	0.125038	0.083648	0.11442	0.142023	0.04802	0.156149	1.10856	0.361184	0.138538	0.163505	0.106313	0.910508	0	0.559447	11.5232	
RY586	0.266104	0.056485	0.127164	0	2.08889	0.73444	0.630813	0.408396	0.369884	0	0	0.225869	0.239794	0.181437	0.136817	0.116315	0.152843	0.070517	0.132131	0.404096	0.2735	0.109494	0.110336	0.103344					

Conventional Compensation workflow



Autospill Compensation Workflow

Choose Selected Detectors

Select multiple detectors by pressing and holding the **shift** key or toggle items by holding **cmd/ctrl** keys. Detectors should at least be equal to the Channels.

- BB515-A
- BB630-A
- BB660-A
- BB700-A
- BB750-A
- BB790-A
- APC-A
- Alexa700-A
- APC-Cy7-A
- BV421-A
- BV480-A
- BV570-A
- BV605-A
- BV650-A
- BV711-A
- BV750-A

Remove Width and Height Detectors.

Choose Selected Detectors Close

Matrix Name: A5 Autospill

Matrix Name: Compensation

Remove Univ Neg Spillover Algorithm Traditional AutoSpill/AutoSpread

Remove the universal negative population (reset to specific negative population)

Parameter	Sample	Comp Name	Negative	Positive
APC-A	CD4	Controls <u>APC</u>	Comp-APC-A	Size/APC-A+
Alexa700-A	CD56 APC R700	Controls <u>Alexa700</u>	Comp-Alexa700-A	Size/Alexa700-A+
BB515-A	CD3 FITC	Controls <u>BB515</u>	Comp-BB515-A	Size/BB515-A+
BB700-A	CD45 PerCP Cy5_5	Controls <u>BB700</u>	Comp-BB700-A	Size/BB700-A+
BUV496-A	CD16	Controls <u>BUV496</u>	Comp-BUV496-A	Size/BUV496-A+
BV480-A	CD8 PB	Controls <u>BV480</u>	Comp-BV480-A	Size/BV480-A+
BV570-A	CD14 BV510	Controls <u>BV570</u>	Comp-BV570-A	Size/BV570-A+
PE-A	CD19	Controls <u>PE</u>	Comp-PE-A	Size/PE-A+

APC-A

Compensation Controls APC Stained Control fcs

SSC-A

FSC-A

APC-A

Spectral Compensation Workflow

Matrix Name: A5 Spectral

Spectral All Detectors Weights Optimize Weights Remove Univ Neg

Spillover Algorithm Traditional AutoSpill/AutoSpread

Confirm gates and control assignments look correct. Double click a graph to edit it.

Parameter	Sample	Comp Name
APC-A	CD4	Comp-APC-A
Alexa700-A	CD56 APC R700	Comp-Alexa700-A
BB515-A	CD3 FITC	Comp-BB515-A
BB700-A	CD45 PerCP Cy5_5	Comp-BB700-A
BUV496-A	CD16	Comp-BUV496-A
BV480-A	CD8 PB	Comp-BV480-A
BV570-A	CD14 BV510	Comp-BV570-A
PE-A	CD19	Comp-PE-A

Spectral All Detectors Weights Optimize Weights Remove Univ Neg

Spillover Algorithm Traditional AutoSpill/AutoSpread

Confirm gates and control assignments look correct. Double click a graph to edit it.

Optimize weights for the all the detectors. This process may take some time to run.


Choose Selected Detectors

Select multiple parameters by pressing and holding the **shift** key or toggle items by holding **cmd/ctrl** keys

- BB515-A
- BB630-A
- BB660-A
- BB700-A
- BB750-A
- BB790-A
- APC-A
- Alexa700-A
- APC-Cy7-A
- BV421-A
- BV480-A
- BV570-A
- BV605-A
- BV650-A
- BV711-A
- BV750-A
- BV786-A

Remove Width and Height Parameters.

Choose Selected Detectors Close



Weights Preview

S#	Detectors	Weights
1	APC-A	1.38981
2	Alexa700-A	2.58485
3	BB515-A	3.87707
4	BB700-A	1.42851
5	BUV496-A	0.6667
6	BV480-A	1.04999
7	BV570-A	0.6667
8	PE-A	1.49993
9	BV421-A	1.13548
10	BUV805-A	1.49993
11	PE Cy5-A	1.49993
12	BV750-A	0.816508
13	BB630-A	2.14286
14	APC-Cy7-A	2.93504
15	BV711-A	0.730027
16	BUV737-A	0.639305
17	BUV615-A	0.46671
18	BB750-A	1.04285
19	BV650-A	0.6667
20	BB790-A	0.6667
21	PECF594-A	1.49993
22	PE-Cy7-A	0.720981
23	BUV563-A	1.49993
24	BV786-A	2.14286
25	BUV396-A	1.0
26	DAPI-A	1.49993
27	PE Cy5.5-A	0.924713
28	BB660-A	2.93504
29	BV605-A	1.32098
30	BUV661-A	0.46671

Import Weights Export Weights Close

Spectral Autospill Compensation Workflow

Matrix Name: A5 Spectral Autospill

Spectral All Detectors Weights Optimize Weights Remove Univ Neg Spillover Algorithm Traditional AutoSpill/AutoSpread

Confirm gates and control assignments look correct. Double click a graph to edit it.

Optimize weights for the all the detectors. This process may take some time to run.

Parameter	Sample	Comp Name
APC-A	CD4	Controls <u>APC</u> Comp-APC-A
Alexa700-A	CD56 APC R700	Controls <u>Alexa700</u> Comp-Alexa700-A
BB515-A	CD3 FITC	Controls <u>BB515</u> Comp-BB515-A
BB700-A	CD45 PerCP Cy5_5	Controls <u>BB700</u> Comp-BB700-A
BUV496-A	CD16	Controls <u>BUV496</u> Comp-BUV496-A
BV480-A	CD8 PB	Controls <u>BV480</u> Comp-BV480-A
BV570-A	CD14 BV510	Controls <u>BV570</u> Comp-BV570-A
PE-A	CD19	Controls <u>PE</u> Comp-PE-A

Choose Selected Detectors

Select multiple parameters by pressing and holding the **shift** key or toggle items by holding **cmd/ctrl** keys



- BB515-A
- BB630-A
- BB660-A
- BB700-A
- BB750-A
- BB790-A
- APC-A
- Alexa700-A
- APC-Cy7-A
- BV421-A
- BV480-A
- BV570-A
- BV605-A
- BV650-A
- BV711-A
- BV750-A
- BV786-A

Remove Width and Height Parameters.

Choose Selected Detectors Close

APC-A

Compensation Controls: APC Stained Control fcs

SSC-A

FSC-A

APC-A

Weights Preview

S#	Detectors	Weights
1	APC-A	1.38981
2	Alexa700-A	2.58485
3	BB515-A	3.87707
4	BB700-A	1.42851
5	BUV496-A	0.6667
6	BV480-A	1.04999
7	BV570-A	0.6667
8	PE-A	1.49993
9	BV421-A	1.13548
10	BUV805-A	1.49993
11	PE Cy5-A	1.49993
12	BV750-A	0.816508
13	BB630-A	2.14266
14	APC-Cy7-A	2.93504
15	BV711-A	0.730027
16	BUV737-A	0.639305
17	BUV615-A	0.46671
18	BB750-A	1.04285
19	BV650-A	0.6667
20	BB790-A	0.6667
21	PECF594-A	1.49993
22	PE-Cy7-A	0.720981
23	BUV563-A	1.49993
24	BV786-A	2.14266
25	BUV396-A	1.0
26	DAPI-A	1.49993
27	PE Cy5.5-A	0.924713
28	BB660-A	2.93504
29	BV605-A	1.32096
30	BUV661-A	0.46671


Import Weights Export Weights Close

Fluorochrome Markers


Fluorochrome	Marker
BUV486	CD16
APC R700	CD56
APC	CD4
PacBlue	CD8
PE	CD19
PerCP-Cy5.5	CD45
FITC	CD3
BV510	CD14


- Stained on lyophilized cells

- Scatter, Doublet and CD45+ gate




A5 Instrument defined







A5- instrument defined




A5- Conventional



A5-Autospill



A5-Spectral



A5- Spectral Autospill

A5 Data

- Identical data set
- Different methods of compensation matrices applied

BD FACSymphony A5 Population Percentages

Parent Population		A5-Acquisition	A5-Traditional	A5-Autospill	A5-Spectral	A5- Spectral autospill	Range
CD45+	CD19+ B cells	11.5%	11.5%	11.5%	11.2%	11.2%	11.2-11.5%
	CD16+	7.54%	7.55%	7.45%	7.53%	7.54%	7.45-7.55%
	Neg	75.4%	75.5%	77.9%	78.3%	77.5%	75.4-78.3%
CD16+	Neutrophils	7.94%	7.94%	7.89%	7.89%	7.80%	7.89-7.94%
	CD56-CD16+ NK Cells	89.7%	89.7%	89.7%	89.7%	89.8%	89.7-89.8%
Neg	T Cells	93.6%	93.7%	93.6%	93.5%	93.8%	93.5-93.8%
	NKT cells	3.24%	3.22%	3.25%	3.13%	3.21%	3.13-3.25%
	CD56+CD16- NK cells	0.75%	0.76%	0.75%	0.77%	0.77%	0.75-0.77%
T cells	CD8 T cells	26.6%	26.6%	25.3%	25.3%	25.5%	25.3-26.6%
	CD4 T cells	69.0%	69.0%	69.2%	69.0%	69.0%	69.0-69.2%

Final thoughts

- All methods of compensation are acceptable
- Conventional and Spectral have a higher chance of bias and human error
- Conventional and Spectral need positive and negative gates defined which could be time consuming
- Autospill could be faster and more efficient for larger panels
- Autospill could take out the bias that conventional has.