


Canto 1 Configuration

APC-Cy7
APC-H7
APC-Alexa 750
APC eFluor 780


APC
Alexa 647
Alexa 633



Alexa 680
Alexa 700
APC-Cy5.5

PerCP-Cy5.5²
685LP 710/50

PerCP²




FITC
GFP
YFP
Alexa488
Venus

PE-Cy5²
635LP 660/20

PE-Cy5.5²
7AAD¹

PE
dTomato
DsRed



PE-CF594
PE-Alexa 610
PI
mCherry

PE-Cy7
PE-Alexa 750

Note: There may be other compatible fluorophores/dyes that are not listed here.
Always check fluorochrome combinations prior to running an experiment!

*Denotes filter change

¹7AAD and PI may be problematic depending on your fluorochrome panel.

²These fluorochromes should not be combined.